



Designation	POWDER FREE NITRILE GLOVES					
Reference(s)		XS	S	M	L	XL
	Blue	3101862	3101863	3101864	3101865	3101866
	White	3101871	3101872	3101873	3101874	-
	Lavender	3101875	3101876	3101877	3101878	-
	Pink	3102752	3102753	3102754	3102755	-
	Blue (200)	3101867	3101868	3101869	3108870	-
Green	3107036	3107037	3107038	3107039		
Visual						
Description	Nitrile examination gloves, powder-free, finger-textured.					
Characteristics and technical datas	<ul style="list-style-type: none"> <li>• Excellent grip.</li> <li>• Allowing the manipulation of silicone.</li> <li>• Conforms to EN 455 &amp; EN 374.</li> </ul>					
Storage	Store protected from the light					
Shelf life	5 years since the date of manufacture					
Indication for use						
Legal information	Medical Device of class I under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking					