



<b>Designation</b>	<b>Brush holders and brushes</b>									
<b>Reference(s)</b>	<table border="1"> <thead> <tr> <th>Designation</th> <th>Reference</th> </tr> </thead> <tbody> <tr> <td>Pack of 10 plastic holders</td> <td>3134142</td> </tr> <tr> <td>Pack of 50 FINES</td> <td>3134138</td> </tr> <tr> <td>Pack of 100 REGULAR</td> <td>3134141</td> </tr> </tbody> </table>		Designation	Reference	Pack of 10 plastic holders	3134142	Pack of 50 FINES	3134138	Pack of 100 REGULAR	3134141
Designation	Reference									
Pack of 10 plastic holders	3134142									
Pack of 50 FINES	3134138									
Pack of 100 REGULAR	3134141									
<b>Product picture</b>										
<b>Product description</b>	<ul style="list-style-type: none"> <li>• Plastic brush holders, cold sterilisation.</li> <li>• Disposable brushes adapting brush holders.</li> </ul>									
<b>Characteristics and composition</b>	<b>Brush holders</b>	<b>Brushes</b>								
	<ul style="list-style-type: none"> <li>- Cold sterilisation</li> <li>- Ergonomic, easy grip</li> <li>- Combines spatula mouth</li> <li>- Length: 147.50 mm (+/- 1 mm)</li> <li>- Polypropylene – 99%</li> <li>- Cadmium free pigments – 1%</li> </ul>	<ul style="list-style-type: none"> <li>- Hard moulded plastic base</li> <li>- Black and soft hair</li> <li>- Polyéthylène – 100%</li> <li>- Black nylon - 100%</li> </ul>								
<b>Packaging</b>	Pack of 10 plastic holders. Bag of 50 or 100 brushes.									
<b>Shelf life</b>	n.a									
<b>Indications of use</b>	<div style="display: flex; align-items: center; gap: 20px;"> <div style="border: 1px solid black; padding: 5px; font-weight: bold;">LOT</div>  </div>									
<b>Legal information</b>	Medical Device of class I CE Marking Manufacturer NOVOBRUSH									