


<b>Designation</b>	<b>Cheek retractors Medibase</b>
<b>Article Reference(s)</b>	Adult: 3135081 Child: 3135080
<b>Product picture</b>	
<b>Product description</b>	<ul style="list-style-type: none"> <li>- Cheek retractors extra resistant and flexible</li> <li>- Easy to set</li> <li>- Ideal for scaling, filling and setting brackets</li> </ul>
<b>Sterilization instruction</b>	<p>Rinse the disinfection materials thoroughly before autoclaving</p> <ul style="list-style-type: none"> <li>- Autoclavable up to 134°C. We recommend the safety program of 121°C</li> <li>- Cycle up to 50x</li> <li>- Please avoid contact with metal instrument when autoclaving</li> </ul>
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>- Made from high quality polycarbonate</li> <li>- Combines high flexibility with very high resistance</li> <li>- Autoclavable up to 121°C during 20 min or 134°C during 18 min</li> <li>- Available in two sizes: adult and child</li> </ul> <p>Store the product in preference in a dry place</p>
<b>Packaging</b>	<ul style="list-style-type: none"> <li>- Box of 2 cheek retractors – adult</li> <li>- Box of 2 cheek retractors – child</li> </ul>
<b>Shelf life</b>	n.a
<b>Indications of use</b>	n.a
<b>Legal information</b>	<p>Medical Device of class I CE Marking Manufacturer: LARIDENT S.R.L.</p>