


Product	Cup brushes							
Reference	<table border="1"> <thead> <tr> <th>Designation</th> <th>Reference</th> </tr> </thead> <tbody> <tr> <td>Nylon</td> <td>3135837</td> </tr> <tr> <td>Silk</td> <td>3135838</td> </tr> </tbody> </table>		Designation	Reference	Nylon	3135837	Silk	3135838
	Designation	Reference						
	Nylon	3135837						
Silk	3135838							
Picture								
Product description	Mini brushes for cleaning and polishing teeth. To be used with contra angle.							
Characteristics & composition	Available in nylon or bristle.							
Packaging	Pack of 50 mini brushes.							
Shelf life	3 years							
Legal informations	Class: IIA. CE marking 0120. SGS UK LIMITED Notified body. Manufacturer: STODDARD & CO.							