

Product	Medibase cup brushes	
Reference	Designation	Reference
	Nylon	3135837
	Silk	3135838
Picture		
Product description	Mini brushes for cleaning and polishing teeth.	
Characteristics	- To be used with a contra-angle: max 10,000 rpm - Single use only	
Packaging	Pack of 50 cup brushes	
Shelf life	5 years	
Indication of use		
Legal informations	Medical Device Class IIA. CE 1639 marking - SGS BELGIUM NV Notified body. Manufacturer: STODDARD & CO.	

Version 2.1 – 10/2025