

<b>Designation</b>	<b>Disposable dental mirrors</b>		
<b>References</b>	<b>Color</b>	White	Assorted (pink, blue, white, green, yellow)
	<b>Reference</b>	3125101	3125100
<b>Product picture</b>			
<b>Product description</b>	Mirrors N ° 5 - Outside diameter: 24 mm Single use only.		
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>- Each mirror is individually wrapped</li> <li>- Each mirror is protected by a plastic "Peel off" to avoid any scratches</li> <li>- Lightweight, with anti-slip system</li> <li>- Mirror: <ul style="list-style-type: none"> <li>• Diameter: 20 mm + 0.15 / 0.20 mm</li> <li>• Thickness: 1.4 mm ± 0.1 mm</li> </ul> </li> </ul>		
<b>Packaging</b>	Economic box of 100 mirrors: white or colors.		
<b>Shelf life</b>	n.a		
<b>Indications of use</b>	n.a		
<b>Legal informations</b>	Class I medical device Manufacturer: CDM CENTER OF EXCELLENCE LTD.		