

Designation	Disposable plastic trays Medibase			
Reference(s)	Size (cm)	With / Without Compartment	Streaked / Smooth	References
	20x10	Without compartment	Streaked	0116227
	18x14			0116228
	28x18			0116229
	20x15		Smooth	3002271
	20x10			3001086
	20x10	With compartment	Streaked	0116204
	18x14			0116205
	28x18			0116203
	28x18		Smooth	3000913
Product picture				
Product description	White single-use non-sterile plastic tray, with compartments and without, smooth and streaked surface. Ideal for setting out dental instruments.			
Characteristics and composition	- Material: High Impact polystyrene (HIPS)			
Packaging	400 trays per carton			
Shelf life	n.a			
Indications of use				
Legal information	n.a			