















Designation	Disposable scalpel															
References	N°10 3131617	N°11 3131618	N°12 3131619	N°15 3131620												
Product picture	 <table border="1" data-bbox="536 855 1374 1120"> <tr> <td>N°10</td> <td></td> <td>10</td> </tr> <tr> <td>N°11</td> <td></td> <td>11</td> </tr> <tr> <td>N°12</td> <td></td> <td>12</td> </tr> <tr> <td>N°15</td> <td></td> <td>15</td> </tr> </table>				N°10		10	N°11		11	N°12		12	N°15		15
N°10		10														
N°11		11														
N°12		12														
N°15		15														
Product description	Disposable scalpel															
Characteristics and composition	<ul style="list-style-type: none"> ✓ Efficiency → Excellent cutting precision ✓ Taken in hand safe and effective → Plastic handle anatomically designed ✓ Resistance → Perfect junction between the blade and the handle ✓ High quality → Stainless Steel blade ✓ Safety → Protection lid ✓ Sterilized by Gamma irradiation ✓ Available in 4 model 															
Packaging	Packaged under individual blister pack															
Shelf life	5 years															
Indications d'utilisation																
Informations légales	Medical device of class IIa Marking CE2460- Certifying body: DNV GL NEMKO PRESAFE AS Manufacturer: RAZORMED INC.															