

Designation	EARLOOP MASKS - PROMASK						
References	White	Blue	Green	Pink	Lavander	Lime	Orange
	3162329	3162328	3162327	3162326	3176259	3176260	3176261
Product picture							
Product description	Medical Masks with elastics triple thickness without latex and without fiberglasses. In compliance with the Typical standard IIR IN 14683: protective masks against the projections of biological liquids.						
Characteristics and composition	▪ Bacterial Filtration Efficiency (BFE): >99%. ▪ Respiratory resistance (ΔP): < 49 Pa/cm². ▪ Splash resistance: > 120 mm Hg. Propreté microbienne (UFC/g) : ≤ 30 ▪ Inner and top layer: Spunbound polypropylene ▪ Filter: Meltblown polypropylene ▪ Dimensions : 17,3 cm x 8,9 cm (± 0,3 cm). ▪ 3 ply						
Packaging	50 masks per box. 10 boxes per carton.						
Shelf life	3 years						
Indications of use							
Legal information	Medical device of class I under MDR 2017/745 Regulation EC Mark – Manufacturer : MEDICOM SAS						

Version 1.2 – 04/08/2025