




Product	DISPOSABLE SYRINGES
Reference	3124778
Picture	
Description	Single-use syringe used for all types of impression materials.
Characteristics & Technical data	<p>Syringe with swivel tip.</p> <ul style="list-style-type: none"> - Volume: 5 ml - Diameter: 12.1 mm ± 0.2 mm - Length: 84 mm ± 1 mm
Storage	Box of 50 syringes
Shell life	 
Legal Information	<p>Class I medical device CE marking Manufacturer : DEPRO AB</p>