



Product	Headrest Covers									
Reference	Colour	White	Light blue	Dark blue	Pink	Lavender	Green	Fuchsia	Lime	Orange
	Ref	141127	141128	141126	3123058	3123060	3123062	3123059	3123061	3123063
Picture										
Description	Single use cellulose headrest covers doubled with polyethylene									
Characteristics & Technical data	<ul style="list-style-type: none"> - Dimensions: 33 x 34 cm - Very durable to tearing paper thanks to its long fibres - The back of the cover is cut halfway to facilitate threading - Cellulose: 26 g/m² - Polyethylene 									
Conditioning	Box of 300 covers (4 dispenser boxes of 75 covers)									
Shelf Life	N.A									
Indications of use										
Legal Information	N.A									