



Designation	Perforated plastic impression tray Medibase										
Reference(s)	Type	Upper			Lower			Upper/Lower			Assorted
	Size	Large	Medium	Small	Large	Medium	Small	Upper left / Lower right	Upper right / Lower Left	Anterior	
	Ref	3135484	3135483	3135490	3135482	3135481	3135489	3135487	3135486	3135488	3135485
	Size A* (mm)	70	65	60	73	70	67	56	56	58	
	Size B* (mm)	60	55	50	60	52	50	50	50	34	
											
Visual											
Description	Disposable plastic impression trays with excellent retention, handling and rigidity.										
Characteristics and composition	<ul style="list-style-type: none"> - Perforated - Excellent retention - Hard plastic - Easy to use - Cold disinfection 										
Packaging	Available in 9 sizes + 1 assorted pack.										
Shelf life	n.a										
Indications of use	<div style="border: 2px solid black; padding: 5px; display: inline-block;"> LOT </div>										
Legal information	Medical device of Class 1 Manufacturer: FOMED Industries Inc										