


Designation	Perforated plastic impression tray Medibase										
Reference(s)	Type	Upper			Lower			Upper/Lower			Assorted
	Size	Large	Medium	Small	Large	Medium	Small	Upper left / Lower right	Upper right/ Lower Left	Anterior	
	Ref	3135484	3135483	3135490	3135482	3135481	3135489	3135487	3135486	3135488	3135485
	Size A* (mm)	70	65	60	73	70	67	56	56	58	
	Size B* (mm)	60	55	50	60	52	50	50	50	34	
											
Visual											
Description	Disposable plastic impression trays with excellent retention, handling and rigidity.										
Characteristics and compositon	<div>- Perforated</div> <div>- Excellent retention</div> <div>- Hard plastic</div> <div>- Easy to use</div> <div>- Cold disinfection</div>										
Packaging	Available in 9 sizes + 1 assorted pack.										
Shelf life	5 years										
Indications of use	<div>LOT</div>										
Legal information	<div>Medical device of Class 1</div> <div>Manufacturer: FOMED Industries Inc</div>										