

Designation	Micro-applicators			
References	SIZE	COLOR	PACKAGING	REFERENCE
	REGULAR	BLUE	Box of 4 tubes each containing 100 micro-applicators	3101006
		BLUE/PINK		3101007
		BLUE/ORANGE		3101008
	FINE	BLACK	Tube of 100 micro-applicators	3101009
SUPER FINE	WHITE	3101010		
Product picture				
Product description	Disposable micro-applicators with foldable head.			
Characteristics and composition	Polypropylene handle, with breakable tip: allows to differentiate the micro-Applicators between them! Nylon bristles. Overall length: 102 mm Tube packaging: hygienic dispensing			
Packaging	REGULAR: box of 4 tubes each containing 100 micro-applicators FINE / SUPER FINE: tubes of 100 micro-applicators			
Shelf life	N/A			
Indication of use	Regular fibers, non-absorbent and lint-free. The solutions are held in suspension between the fibers to prevent dripping and allow perfect application. Easy access to all areas of application thanks to the foldable head. Disposable: discard after use. For use with silanes, adhesives, etching gel, hemostatic solutions, cavity bottoms and plaque developers. Store away from light and dust in a cool, dry place.			
Legal informations	Class I medical device Manufacturer: AKZENTA INTERNATIONAL CHINA			