

Product	Mixing-tips and intra-oral tips								
Reference, Visual and Dimensions		For light addition silicone		For bite registration	For bonding and cement composite		For automatic mixiing	For composites, crowns and briges	For etching
									
	Reference s	3149958	3150588	3149957	3149960	3149956	3180067	3149959	3180066
	Packaging (pieces)	100		100	50		50	50	50
	Ratio	1:1	1:1	1:1	1:1	1:1	5:1	10 : 1	NA
	Base length	1,8 cm	NA	1,9 cm	1,9 cm	1,8 cm	1,2 cm	2 cm	1,1 cm
	Total length	7 cm	2.5 cm	8,9 cm	6,3 cm	3,4 cm	5 cm	7,5 cm	3 cm
Description	Complete range of single use mixing tips and intra oral tips.								
Packaging	50 or 100 tips per bags.								
Indications of use									
Legal information	Class I medical device in accordance with EU Regulation 2017/745 (MDR) CE marked Manufacturer : - KUNSHAN NSJ MEDICAL EQUIPMENT CO, LTD. (ref 3149958, 3149957, 3149960, 3149956, 3180067, 3149959, 3180066) - SEILGLOBAL CO. LTD. (ref 3150588)								

Version 5.0 – 11/2025