



Designation	Non woven swabs, non sterile Medibase											
Reference(s)	3131819											
Product picture												
Product description	<ul style="list-style-type: none"> • Do not realise fibers when cleaning wounds • Good skin tolerance, even for sensitive skin • Compatible with antiseptic solutions and ether. • Homogeneous structure allowing a sensation of softness in contact with the skin and a low risk of fraying. • Absence of allergenic chemical binder <table border="1"> <thead> <tr> <th>Weight/m²</th> <th>Number of thickness</th> <th>Folded dimensions</th> <th>Unfolded dimensions</th> </tr> </thead> <tbody> <tr> <td>40 g</td> <td>4</td> <td>5 x 5 cm</td> <td>10 x 10 cm</td> </tr> </tbody> </table>				Weight/m ²	Number of thickness	Folded dimensions	Unfolded dimensions	40 g	4	5 x 5 cm	10 x 10 cm
Weight/m ²	Number of thickness	Folded dimensions	Unfolded dimensions									
40 g	4	5 x 5 cm	10 x 10 cm									
Characteristics and composition	Composition: 70% viscose - 30% polyester Possible to sterilize with: Vapour, Ethylene Oxide, Radiation											
Packaging	Box of 1 000 swabs (10 packs of 100 swabs)											
Shelf life	Not applicable											
Indications of use												
Legal information	Medical Device of class I CE Marking Manufacturer: SYLAMED											