


Product	PHOTO RETRACTORS
Reference	Adult : 3135082 Infant : 3135084
Picture	
Description	With an extra-long handle, it is ideal for clinical photography
Sterilization advice	Rinse disinfectants thoroughly before autoclaving. - Autoclavable at 134°C for 18 min. - Cycle: up to 50 x - Do not put metal in contact in the autoclave
Characteristics & Technical data	- Non-reflective - Made of high-quality polycarbonate - Sterilizable in autoclave at 134 ° C for 18 min - Available in two sizes: adult and child Store the product preferably in a dry place.
Storage	- Box of 2 adult cheeks retractors - Box of 2 child cheeks retractors
Garantee	More or less 5 years depending on the number of sterilization cycles.
Shell life	n.a
Legal Information	Class I medical CE marking Manufacturer: LARIDENT S.R.L.