


Designation	POWDER FREE VINYL GLOVES MEBIBASE																									
Reference	Size	XS	S	M	L	XL																				
	Ref	3145827	3145828	3145829	3145830	3145831																				
Visual	 																									
Description	Disposable powder free vinyl gloves																									
Characteristics and composition	- Ambidextrous - Powder-free: reduced risk of allergy - High resistance - Colour: transparent - Rolled cuffs prevent gloves from splitting tearing.																									
	<table><tr><th>Properties</th><th>Size</th><th>In-House standard</th></tr><tr><td>Length (mm)</td><td>All sizes</td><td>Min. 240</td></tr><tr><td rowspan="5">Palm width (mm)</td><td>XS</td><td>75 +/- 5</td></tr><tr><td>S</td><td>85 +/- 5</td></tr><tr><td>M</td><td>95 +/- 5</td></tr><tr><td>L</td><td>105 +/- 5</td></tr><tr><td>XL</td><td>115 +/- 5</td></tr><tr><td>Thickness (mm)</td><td>All sizes</td><td>Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03</td></tr></table>						Properties	Size	In-House standard	Length (mm)	All sizes	Min. 240	Palm width (mm)	XS	75 +/- 5	S	85 +/- 5	M	95 +/- 5	L	105 +/- 5	XL	115 +/- 5	Thickness (mm)	All sizes	Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03
	Properties	Size	In-House standard																							
	Length (mm)	All sizes	Min. 240																							
	Palm width (mm)	XS	75 +/- 5																							
S		85 +/- 5																								
M		95 +/- 5																								
L		105 +/- 5																								
XL		115 +/- 5																								
Thickness (mm)	All sizes	Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03																								
<table><tr><th>Characteristics</th><th>In-House standard</th></tr><tr><td>Weight (g)</td><td>4.5±0.3g (Size M)</td></tr><tr><td>Powder Level</td><td>≤2mg/glove</td></tr><tr><td>Tensile Strength(N) Before & after aging</td><td>3.6</td></tr></table>						Characteristics	In-House standard	Weight (g)	4.5±0.3g (Size M)	Powder Level	≤2mg/glove	Tensile Strength(N) Before & after aging	3.6													
Characteristics	In-House standard																									
Weight (g)	4.5±0.3g (Size M)																									
Powder Level	≤2mg/glove																									
Tensile Strength(N) Before & after aging	3.6																									

	Complies with the standards:	EN 455-1:2020; EN 455-2:2015; EN 455-3:2015; EN 455-4:2009; ISO 10993-1:2018; ISO 10993-10:2010; ISO 10993-5:2009; EN ISO 21420:2020; EN ISO 374-1:2016+A1:2018; EN ISO 374-4:2019; EN ISO 374-5:2016;
Packaging	Box containing 100 gloves. 10 boxes per carton.	
Shelf life	5 years from the manufacturing date	
Indications of use		
Legal information	Medical device of class I under MDR 2017/745 regulation + EPI CAT III CE 2777 Marking Notified body : SATRA Technology Europe Ltd. Manufacturer: ANHUI INTCO MEDICAL PRODUCTS CO.,Ltd.	

Version 1.0

Updated 10/12/2021 by MF