


Designation	POWDER FREE VINYL GLOVES MEBIBASE																									
Reference	<table><tr><td>Size</td><td>XS</td><td>S</td><td>M</td><td>L</td><td>XL</td></tr><tr><td>Ref</td><td>3145827</td><td>3145828</td><td>3145829</td><td>3145830</td><td>3145831</td></tr></table>						Size	XS	S	M	L	XL	Ref	3145827	3145828	3145829	3145830	3145831								
Size	XS	S	M	L	XL																					
Ref	3145827	3145828	3145829	3145830	3145831																					
Visual																										
Description	Powder-free, single-use vinyl (Polyvinyl Chloride) examination gloves																									
Characteristics and composition	- Ambidextrous - Powder-free: reduced risk of allergy - High resistance - Colour: transparent - Rolled cuffs prevent gloves from splitting tearing.																									
	<table><tr><th>Properties</th><th>Size</th><th>In-House standard</th></tr><tr><td>Length (mm)</td><td>All sizes</td><td>Min. 240</td></tr><tr><td rowspan="5">Palm width (mm)</td><td>XS</td><td>75 +/- 5</td></tr><tr><td>S</td><td>85 +/- 5</td></tr><tr><td>M</td><td>95 +/- 5</td></tr><tr><td>L</td><td>105 +/- 5</td></tr><tr><td>XL</td><td>115 +/- 5</td></tr><tr><td>Thickness (mm)</td><td>All sizes</td><td>Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03</td></tr></table>						Properties	Size	In-House standard	Length (mm)	All sizes	Min. 240	Palm width (mm)	XS	75 +/- 5	S	85 +/- 5	M	95 +/- 5	L	105 +/- 5	XL	115 +/- 5	Thickness (mm)	All sizes	Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03
	Properties	Size	In-House standard																							
	Length (mm)	All sizes	Min. 240																							
Palm width (mm)	XS	75 +/- 5																								
	S	85 +/- 5																								
	M	95 +/- 5																								
	L	105 +/- 5																								
	XL	115 +/- 5																								
Thickness (mm)	All sizes	Finger: 0.08 +/- 0.03 Palm: 0.07 +/- 0.03 Cuff: 0.06 +/- 0.03																								
Nature of the additives, curing agents or accelerators used :		DOTP, Calcium-Zinc Complex Stabilizer, Polyurethane																								
Quantity of additives, curing agents or accelerators used :		DOTP(40%-43%) Calcium-Zinc Complex Stabilizer (0.4%-0.6%) Polyurethane (2%-3%)																								

	Characteristics	In-House standard
	Weight (g)	4.5±0.3g (Size M)
	Powder Level	≤2mg/glove
	Tensile Strength(N) Before & after aging	3.6
	Complies with the standards:	EN 455-1:2020; EN 455-2:2015; EN 455-3:2015; EN 455-4:2009; ISO 10993-1:2018; ISO 10993-10:2010; ISO 10993-5:2009; EN ISO 21420:2020; EN ISO 374-1:2016+A1:2018; EN ISO 374-4:2019; EN ISO 374-5:2016;
Packaging	Box containing 100 gloves. 10 boxes per carton.	
Shelf life	5 years from the manufacturing date	
Indications of use		
Legal information	Medical device of class I under MDR 2017/745 regulation + EPI CAT III CE 2777 Marking Notified body: SATRA Technology Europe Ltd. Manufacturer: ANHUI INTCO MEDICAL PRODUCTS CO.,Ltd.	

Version 2.0

Updated 13/12/2022 by RG