



<b>Designation</b>	<b>POWDER FREE LATEX EXAMINATION GLOVES SENSI +</b>																
<b>Reference</b>	XS – 5/6	S – 6/7	M – 7/8	L – 8/9	XL – 9/10												
	3135154	3135155	3135156	3135157	3135158												
<b>Visual product</b>																	
<b>Product description</b>	Disposable powder free, micro-textured, natural latex rubber examination gloves. Ambidextrous, single use.																
<b>Characteristics and composition</b>	<table border="1"> <thead> <tr> <th></th> <th>XS-5/6</th> <th>S-6/7</th> <th>M-7/8</th> <th>L-8/9</th> <th>XL-9/10</th> </tr> </thead> <tbody> <tr> <td>Width of the palm (mm)</td> <td>70-79</td> <td>80-89</td> <td>90-99</td> <td>100-109</td> <td>110-119</td> </tr> </tbody> </table>						XS-5/6	S-6/7	M-7/8	L-8/9	XL-9/10	Width of the palm (mm)	70-79	80-89	90-99	100-109	110-119
		XS-5/6	S-6/7	M-7/8	L-8/9	XL-9/10											
Width of the palm (mm)	70-79	80-89	90-99	100-109	110-119												
<ul style="list-style-type: none"> <li>– Excellent resistance and elasticity.</li> <li>– Designed for excellent fit and comfort.</li> <li>– Maximum dexterity and tactile sensitivity (even with wet and slippery materials).</li> <li>– Rolled cuffs prevent gloves from splitting tearing.</li> <li>– Comply with standards EN 455-1/2/3/4, EN 420:2003+A1:2009, EN 374-1:2016 EN374-5:2016</li> </ul>																	
<b>Packaging</b>	Box included 100 gloves ( $\pm$ 2 gloves). 10 boxes per carton.																
<b>Shelf life</b>	5 years from the date of manufacture.																
<b>Indication for use</b>																	
<b>Legal information</b>	Medical Device of class I under MDR (EU) 2017/745 + + EPI Cat. III CE Marking 0598. Notified body: SGS FIMKO OY. Manufacturer: MAXTER GLOVE MANUFACTURING Sdh Bhd.																

Version 1.1  
Update 2022/03/04 by EH