



| Product Designation | DISPOSABLE SYRINGES FOR SINGLE-USE | | | | | | | | | | | | | | | | | | |
|--|--|----------------------------|---------|---------|---------|------------|----------------------------|--------|--------|-----|------|--------|-----|-------|--------|-----|-------|------|----|
| Reference | Content | 2,5 ML | 5 ML | 10 ML | 20 ML | | | | | | | | | | | | | | |
| | Reference | 3125103 | 3125105 | 3125102 | 3125104 | | | | | | | | | | | | | | |
| Product picture |  | | | | | | | | | | | | | | | | | | |
| Product description | <ul style="list-style-type: none"> ▪ Disposable syringes for single-use for endodontic canal irrigation | | | | | | | | | | | | | | | | | | |
| Characteristics and composition | <ul style="list-style-type: none"> ▪ Luer-lock thread for a safe fixation of the needle ▪ Latex free piston ▪ Sterile individual package | | | | | | | | | | | | | | | | | | |
| Packaging | <table border="1"> <thead> <tr> <th>Content</th> <th>Graduation</th> <th>Number of syringes per box</th> </tr> </thead> <tbody> <tr> <td>2.5 ML</td> <td>0.1 ML</td> <td>100</td> </tr> <tr> <td>5 ML</td> <td>0.2 ML</td> <td>100</td> </tr> <tr> <td>10 ML</td> <td>0.5 ML</td> <td>100</td> </tr> <tr> <td>20 ML</td> <td>1 ML</td> <td>80</td> </tr> </tbody> </table> | | | | Content | Graduation | Number of syringes per box | 2.5 ML | 0.1 ML | 100 | 5 ML | 0.2 ML | 100 | 10 ML | 0.5 ML | 100 | 20 ML | 1 ML | 80 |
| Content | Graduation | Number of syringes per box | | | | | | | | | | | | | | | | | |
| 2.5 ML | 0.1 ML | 100 | | | | | | | | | | | | | | | | | |
| 5 ML | 0.2 ML | 100 | | | | | | | | | | | | | | | | | |
| 10 ML | 0.5 ML | 100 | | | | | | | | | | | | | | | | | |
| 20 ML | 1 ML | 80 | | | | | | | | | | | | | | | | | |
| Shelf life | 5 years from the manufactured date | | | | | | | | | | | | | | | | | | |
| Indications of use |  | | | | | | | | | | | | | | | | | | |
| Legal information | Medical Device of class IIA CE Marking- 0197-TUV RHEINLAND Manufacturer: CDM CENTER OF EXCELLENCE LTD | | | | | | | | | | | | | | | | | | |