



<b>Product Designation</b>	<b>DISPOSABLE SYRINGES FOR SINGLE-USE</b>																		
<b>Reference</b>	<table border="1"> <tr> <td><b>Content</b></td> <td>2,5 ML</td> <td>5 ML</td> <td>10 ML</td> <td>20 ML</td> </tr> <tr> <td><b>Reference</b></td> <td>3125103</td> <td>3125105</td> <td>3125102</td> <td>3125104</td> </tr> </table>	<b>Content</b>	2,5 ML	5 ML	10 ML	20 ML	<b>Reference</b>	3125103	3125105	3125102	3125104								
<b>Content</b>	2,5 ML	5 ML	10 ML	20 ML															
<b>Reference</b>	3125103	3125105	3125102	3125104															
<b>Product picture</b>																			
<b>Product description</b>	<ul style="list-style-type: none"> <li>Disposable syringes for single-use for endodontic canal irrigation</li> </ul>																		
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>Luer-lock thread for a safe fixation of the needle</li> <li>Latex free piston</li> <li>Sterile individual package</li> </ul>																		
<b>Packaging</b>	<table border="1"> <thead> <tr> <th>Content</th> <th>Graduation</th> <th>Number of syringes per box</th> </tr> </thead> <tbody> <tr> <td>2.5 ML</td> <td>0.1 ML</td> <td>100</td> </tr> <tr> <td>5 ML</td> <td>0.2 ML</td> <td>100</td> </tr> <tr> <td>10 ML</td> <td>0.5 ML</td> <td>100</td> </tr> <tr> <td>20 ML</td> <td>1 ML</td> <td>80</td> </tr> </tbody> </table>				Content	Graduation	Number of syringes per box	2.5 ML	0.1 ML	100	5 ML	0.2 ML	100	10 ML	0.5 ML	100	20 ML	1 ML	80
Content	Graduation	Number of syringes per box																	
2.5 ML	0.1 ML	100																	
5 ML	0.2 ML	100																	
10 ML	0.5 ML	100																	
20 ML	1 ML	80																	
<b>Shelf life</b>	5 years from the manufactured date																		
<b>Indications of use</b>																			
<b>Legal information</b>	Medical Device of class IIA CE Marking- 0197-TUV RHEINLAND Manufacturer: CDM CENTER OF EXCELLENCE LTD																		