

## **TECHNICAL SHEET**

Designation	Hydrophilic gauze swabs sterile Medibase
Reference(s)	3135500
Product picture	A COC. Statements of grave and statements of the control of the co
Product description	<ul> <li>Dimensions folded: 5*5 cm</li> <li>Number of threads: 17</li> <li>Number of plies: 12</li> <li>Sterilisation: Ethylene Oxide</li> </ul>
Characteristics and compositon	<ul> <li>Dimensions folded: 5*5 cm</li> <li>Number of threads: 17</li> <li>Number of plies: 12</li> <li>Sterilisation: Ethylene Oxide</li> </ul>
Packaging	Carton of 100 swabs (50 sterile pouches of 2 swabs)
Shelf life	5 years
Indications of use	STERILEEO
Legal information  Version 1.1	Medical Device of class I CE 0123 Marking- Certifying body: TUV PRODUCT SERVICE Manufacturer: FOMED Industries Inc Updated 2021/12/08 by Maria FRENZ

Version 1.1

Updated 2021/12/08 by Maria FRENZ