


Designation	Sterile needles Medibase (100)					
Reference(s)		Description	G	Ø	Longueur	Reference
	Périapicales	Extra-fine	30 G	0,30	16 mm	3114384
		Extra-fine	30 G	0,30	21 mm	3114385
		Fine	27 G	0,40	8 mm	3114386
			25G	0,50	36 mm	3114387
			27 G	0,40	21 mm	3114388
	Intraseptales	Fine	27 G	0,40	16 mm	3114389
	Intraligamentaire	Extra-fine	30 G	0,30	12 mm	3114390
	Tronculaires	Extra-fine	30 G	0,30	25 mm	3114391
Visual						
Description	Sterile needles for single use					
Features and Technical Data	- Excellent quality, with silicone coating - Triple bevel for a gentle injection - Sterile with sterilization control ring - With bevel indicator - Sterilized with ethylene oxide - Cannula: stainless steel - Plastic part: polypropylene (PP) - Lubricant: silicone					
Packaging	Box of 100 needles					
Lifetime	5 years from the date of manufacture					
Indications of use						
legal informations	Class IIA medical device CE Marking 0426. ITALCERT Notify body Manufacturer: TECNOFAR SPA					