



| | | | | |
|--|---|-----------------------|-----------------------|-----------------------|
| Designation | Sterile non-woven swabs Medibase | | | |
| Reference(s) | Size | 5x5 cm | 5x5 cm | 7,5x7,5 cm |
| | Box of | 150 swabs | 250 swabs | 80 swabs |
| | Internal conditioning | 75 pouches of 2 swabs | 50 pouches of 5 swabs | 40 pouches of 2 swabs |
| | Reference | 3135502 | 3135501 | 3135503 |
| Product picture |  | | | |
| Product description | Non-woven swabs Medibase packed under sterile pouches | | | |
| Characteristics and composition | Swabs 60% viscose / 40% polyester <ul style="list-style-type: none"> • Number of plies : 4 • Sterilisation : ethylene oxide | | | |
| Packaging | Available in two sizes (5x5 cm and 7,5x7,5 cm) and three conditioning | | | |
| Shelf Life | 5 years | | | |
| Indications of use |  | | | |
| Legal information | Medical Office of class I CE 0123 Marking. Certifying body: TUV PRODUCT SERVICE Manufacturer: FOMED Industries Inc | | | |