


Designation	Sterilisation pouches Medibase									
Reference(s)	SELF SEALING POUCHES									
	Size	8,5 x 16,5cm	7 x 25,5cm	9 x 25,5cm	14 x 28cm	25,5 x 38cm	5,7 x 13cm	18 x 33cm	9 x 13cm	30 x 45 cm
	Ref	3135469	3135470	3135471	3135472	3135473	3135474	3135475	3135476	3135477
	HEAT SEALING POUCHES									
	Size	7 x 25,5cm			9 x 25,5cm			20 x 33cm AVEC SOUFFLET DE 5,5cm		
	Ref	3135479			3135480			3135478		
Product picture										
Product description	Sterilization bags recommended for autoclaving or chemiclave.									
Characteristics and composition	Paper	The basis paper weight is $\pm 60\text{g/m}^2$, fitting the requirements of sterile packaging manufactured paper which list on ISO11607-1, 2 and EN868-5; According to the pouch designs, the printed inscription on the rolls is navy blue (Pantone298U), printing steam and EO indicator.								
	Film	Lucid wathet blue hyperthermal PET/CPP complex film, thickness is $\pm 0,052\text{mm}$, weight is 56 g/m^2 , conforming to the related requirements of EN868-5.								
	Ink	Using ink; the min printing area of indicator should be 100mm^2 , may not fade or become stained.								

Packaging	Box of 200 pouches, available in 10 sizes
Shelf life	5 years from manufacturing date
Indications of use	
Legal information	Medical Device of class I CE Marking Manufacturer: FOMED Industries Inc

Version 2.2

Updated 2022/01/03 by Maria FRENZ