


Designation	STERILIZATION ROLLS						
References	5 cm	7,5 cm	10 cm	15 cm	20 cm	25 cm	20 cm wide-sided depth 5,5 cm
	3134962	3134963	3134964	3134965	3134966	3134967	3134968
Visual							
Description	Sterilization pouches in roll of 200 m, available in 6 different widths + one size wide-sided (roll of 100m).						
Characteristics	<p>Medical paper: 60 grams PET films: 12/30 of blue color (12/40 for the girdle with bellows)</p> <p>Every roll is individually packed under plastic film.</p> <ul style="list-style-type: none"> ▪ Temperature of sealing of the bag: 170°C-200°C ▪ Steam sterilization <ul style="list-style-type: none"> ○ 121°C in 134°C ○ Indicators turn pink to brown ▪ Sterilization with oxide of Ethylene <ul style="list-style-type: none"> ○ 30°C in 54°C ○ Indicators turn green to yellow ▪ In accordance with the standards: <ul style="list-style-type: none"> ○ 868:5 ○ ISO 11607: 1-2 ○ ISO 1140-1 						
Storage	KEEP IN DRY PLACE, KEEP AWAY FROM SUNLIGHT DO NOT USE IF PACKAGING IS DAMAGED						
Shelf-life	5 years from the date of manufacture						
Legal information	Medical device: class I CE MARKING Manufacturer: E-LINE						