

|                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Designation</b>                      | <b>TIE-ON MASKS - PROMASK</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reference(s)</b>                     | Blue<br>3162330                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product picture</b>                  |                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product description</b>              | Medical tie-on masks, 3 plies, latex and fiberglass free. Complies with EN 14683 norm Type IIR: protective masks <b>against the projections of biological liquids.</b>                                                                                                                                                                                                                                                                                                              |
| <b>Characteristics and compositions</b> | <ul style="list-style-type: none"> <li>▪ Bacterial Filtration Efficiency (BFE): &gt;99%.</li> <li>▪ Respiratory resistance (<math>\Delta P</math>): &lt; 49 Pa/cm<sup>2</sup>.</li> <li>▪ Splash resistance: &gt; 120 mm Hg. Propreté microbienne (UFC/g) : <math>\leq 30</math></li> <li>▪ Inner and top layer: Spunbound polypropylene</li> <li>▪ Filter: Meltblown polypropylene</li> <li>▪ Dimensions : 17,3 cm x 8,9 cm (<math>\pm 0,3</math> cm).</li> <li>▪ 3 ply</li> </ul> |
| <b>Packaging</b>                        | 50 masks per box.<br>6 boxes per carton.                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Shelf life</b>                       | 3 years                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Indications of use</b>               |                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Legal information</b>                | Medical Device of class I under MDR 2017/745 regulation<br>EC Mark - Manufacturer MEDICOM SAS                                                                                                                                                                                                                                                                                                                                                                                       |

Version 1.1 – 14/08/2025