



Designation	Mixing spatulas Medibase
Reference(s)	3021027
Product picture	
Product description	Double-ended mixing spatula in Acrylonitrile butadiene styrene (ABS)
Characteristics and composition	<ul style="list-style-type: none"> - Aspect : Cylindrical hard plastic flat tip - Color : White - Length : 11.74 cm - Diameter : 0.41 cm - MEDIBASE logo embossed on each spatula
Packaging	Lot of 100 spatulas in plastic bag
Shelf life	N/A
Indications of use	
Legal information	N/A

Version 2.1

Updated 2016/06/27 by Scarlet Ceausoiu