



Designation	EARLOOP MASKS				
Reference(s)	Blue	White	Green	Pink	Orange
	0116756	0124660	0116757	3010233	3017821
Product picture					
Product description	<ul style="list-style-type: none"> <li>Medical Masks with elastics triple thickness without latex and without fibreglasses. In compliance with the Typical standard IIR IN 14683: protective masks against the projections of biological liquids.</li> </ul>				
Characteristics and composition	<ul style="list-style-type: none"> <li>Bacterial Filtration Efficiency (BFE): &gt;99%.</li> <li>Respiratory resistance (<math>\Delta P</math>): &lt; 49 Pa/cm<sup>2</sup>.</li> <li>Splash resistance: &gt; 120 mm Hg.</li> <li>Integrated nose piece.</li> <li>Dimensions: 8.9 cm x 17.3 cm.</li> <li>Dimensions : 17,3 cm x 8,9 cm (<math>\pm 0,3</math> cm).</li> <li>3 ply of 15 mm (+/- 3 mm).</li> </ul>				
Packaging	Box of 50 masks. 10 boxes per carton.				
Shelf life	NA				
Indications of use					
Legal information	Medical device of class I. Manufacturer: CFPM.				

Version 2.0  
Updated 2016/05/13 by SMA