

Designation	POWDER FREE NITRILE GLOVES					
Reference(s)		XS	S	M	L	XL
	Blue (200)	3135177	3135178	3135179	3135180	-
	Blue (100)	3135172	3135173	3135174	3135175	3135176
	White (100)	3135181	3135182	3135183	3135184	-
	Lavander (100)	3135185	3135186	3135187	3135188	-
	Pink (100)	3135189	3135190	3135191	3135192	-
	Green (100)	3135193	3135194	3135195	3135196	-
Visual						
Description	Nitrile examination gloves, powder-free, finger-textured.					
Characteristics and technical datas	• Excellent grip. • Allowing the manipulation of silicone. • Conforms to EN 455: 1-4 & EN ISO 374-1;4;5.					
Storage	Store protected from the light					
Shelf life	5 years since the date of manufacture					
Indication for use	EN ISO 374-5:2016  VIRUS EN ISO 374-1:2016 Type B  JKT        40°C 104°F max 					
Legal information	Medical Device of class I under MDR directive (EU) 2017/745 PPE Cat III CE Marking 0598					