

Designation	POWDER FREE NITRILE GLOVES GREEN																
Reference(s)	Size	XS	S	M	L												
	Reference	3135193	3135194	3135195	3135196												
	UDI-DI	3665318009472	3665318009465	3665318009458	3665318009441												
Visual																	
Description	Nitrile examination gloves, powder-free. Ambidextrous, finger-textured. Non sterile, beaded cuff																
Characteristics and technical datas	• Excellent grip. • Allowing the manipulation of silicone. • Conforms to EN 455: 1-4 & EN ISO 374-1;4;5. • Food contact approved																
	<table><tr><th>Size</th><th>Glove weight (g)</th><th>Length (mm)</th></tr><tr><td>XS</td><td>2.9 ± 0.2</td><td rowspan="4">Min. 240</td></tr><tr><td>S</td><td>3.2 ± 0.2</td></tr><tr><td>M</td><td>3.4 ± 0.2</td></tr><tr><td>L</td><td>3.7 ± 0.2</td></tr></table>		Size	Glove weight (g)	Length (mm)	XS	2.9 ± 0.2	Min. 240	S	3.2 ± 0.2	M	3.4 ± 0.2	L	3.7 ± 0.2			
	Size	Glove weight (g)	Length (mm)														
	XS	2.9 ± 0.2	Min. 240														
	S	3.2 ± 0.2															
	M	3.4 ± 0.2															
	L	3.7 ± 0.2															
	<table><tr><th rowspan="3">Thickness (mm)</th><td>Finger</td><td>0.10</td></tr><tr><td>Palm</td><td>0.05</td></tr><tr><td>Cuff</td><td>0.05</td></tr></table>		Thickness (mm)	Finger	0.10	Palm	0.05	Cuff	0.05								
	Thickness (mm)	Finger		0.10													
		Palm		0.05													
Cuff		0.05															
<table><tr><th></th><th>Before Aging</th><th>After aging</th></tr><tr><td>Tensile strength (MPa)</td><td>Min. 14</td><td>Min. 14</td></tr><tr><td>Elongation (%)</td><td>Min. 500</td><td>Min. 400</td></tr><tr><td>Force at break (N)</td><td>Min. 6</td><td>Min. 6</td></tr></table>				Before Aging	After aging	Tensile strength (MPa)	Min. 14	Min. 14	Elongation (%)	Min. 500	Min. 400	Force at break (N)	Min. 6	Min. 6			
	Before Aging	After aging															
Tensile strength (MPa)	Min. 14	Min. 14															
Elongation (%)	Min. 500	Min. 400															
Force at break (N)	Min. 6	Min. 6															
<table><tr><td>Protein Level</td><td>Not applicable</td></tr><tr><td>Powder Level</td><td>Max.2 mg/gant</td></tr></table>			Protein Level	Not applicable	Powder Level	Max.2 mg/gant											
Protein Level	Not applicable																
Powder Level	Max.2 mg/gant																
Packaging	Box of 100 gloves Carton of 10 boxes																
Storage	Store protected from the light, water, humidity, and below 40°C																

Shelf life	5 years since the date of manufacture
Indication for use	EN ISO 374-5:2016  VIRUS EN ISO 374-1:2016 Type B  JKT    NON STERILE     40°C 
Legal information	Medical Device of class I under MDR REGULATION (EU) 2017/745 PPE Cat III UNDER PPE REGULATION (EU) 2016/425 - CE Marking 0598

Version 2.0 created by SMA

Date : 12/01/2023