





Designation	POWDER FREE NITRILE GLOVES GREEN																
Reference(s)	Size	XS	S	M	L												
	Reference	3135189	3135190	3135191	3135192												
	UDI-DI	3665318013592	3665318013585	3665318013578	3665318013547												
Visual																	
Description	Nitrile examination gloves, powder-free. Ambidextrous, finger-textured. Non sterile, beaded cuff																
Characteristics and technical datas	<ul style="list-style-type: none"> • Excellent grip. • Allowing the manipulation of silicone. • Conforms to EN 455: 1-4 & EN ISO 374-1;4;5. • Food contact approved 																
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Size</th> <th>Glove weight (g)</th> <th>Length (mm)</th> </tr> </thead> <tbody> <tr> <td>XS</td> <td>2.9 ± 0.2</td> <td rowspan="4" style="text-align: center; vertical-align: middle;">Min. 240</td> </tr> <tr> <td>S</td> <td>3.2 ± 0.2</td> </tr> <tr> <td>M</td> <td>3.4 ± 0.2</td> </tr> <tr> <td>L</td> <td>3.7 ± 0.2</td> </tr> </tbody> </table>					Size	Glove weight (g)	Length (mm)	XS	2.9 ± 0.2	Min. 240	S	3.2 ± 0.2	M	3.4 ± 0.2	L	3.7 ± 0.2
	Size	Glove weight (g)	Length (mm)														
	XS	2.9 ± 0.2	Min. 240														
	S	3.2 ± 0.2															
M	3.4 ± 0.2																
L	3.7 ± 0.2																
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="3">Thickness (mm)</th> <th>Finger</th> <th>0.10</th> </tr> </thead> <tbody> <tr> <th>Palm</th> <th>0.05</th> </tr> <tr> <th>Cuff</th> <th>0.05</th> </tr> </tbody> </table>					Thickness (mm)	Finger	0.10	Palm	0.05	Cuff	0.05						
Thickness (mm)	Finger	0.10															
	Palm	0.05															
	Cuff	0.05															
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th>Before Aging</th> <th>After aging</th> </tr> </thead> <tbody> <tr> <td>Tensile strength (MPa)</td> <td>Min. 14</td> <td>Min. 14</td> </tr> <tr> <td>Elongation (%)</td> <td>Min. 500</td> <td>Min. 400</td> </tr> <tr> <td>Force at break (N)</td> <td>Min. 6</td> <td>Min. 6</td> </tr> </tbody> </table>						Before Aging	After aging	Tensile strength (MPa)	Min. 14	Min. 14	Elongation (%)	Min. 500	Min. 400	Force at break (N)	Min. 6	Min. 6	
	Before Aging	After aging															
Tensile strength (MPa)	Min. 14	Min. 14															
Elongation (%)	Min. 500	Min. 400															
Force at break (N)	Min. 6	Min. 6															
<table border="1" style="width: 100%; border-collapse: collapse;"> <tbody> <tr> <td>Protein Level</td> <td>Not applicable</td> </tr> <tr> <td>Powder Level</td> <td>Max.2 mg/gant</td> </tr> </tbody> </table>					Protein Level	Not applicable	Powder Level	Max.2 mg/gant									
Protein Level	Not applicable																
Powder Level	Max.2 mg/gant																
Packaging	Box of 100 gloves Carton of 10 boxes																
Storage	Store protected from the light, water, humidity, and below 40°C																

Shelf life	5 years since the date of manufacture
Indication for use	<p>EN ISO 374-5:2016 EN ISO 374-1:2016 Type B</p>  <p>VIRUS JKT LATEX NON STERILE</p>  <p>40°C</p>  <p>AQL 1.5</p>
Legal information	<p>Medical Device of class I under MDR REGULATION (EU) 2017/745</p> <p>PPE Cat III UNDER PPE REGULATION (EU) 2016/425 - CE Marking 0598</p>

Version 2.0 created by SMA

Date : 12/01/2023