


Designation	Saliva ejectors – fixed tips		
References	Colours	transparent	white
	Reference	3135200	3135201
Product Picture			
Product description	<p>Easy shaping because of copper wire reinforcement</p> <ul style="list-style-type: none"> - Tips : small, smooth and rounded - Ideal for all uses, even on irritated walls - Diameter : 6,5 mm - Length : 15 cm 		
Characteristics and composition	Made from PVC		
Packaging	<ul style="list-style-type: none"> - Store at room temperature in a dry and ventilated area - Plastic bag of 100 units 		
Shelf life	5 years from the date of manufacture		
Indications of use			
Legal information	<p>Medical Device of class IIA CE marking- 0051-IMQ S.p.A Manufacturer: ASA DENTAL</p>		

Designation	Saliva ejectors – removable tips					
Reference(s)	Colours	Translucent	Blue	Green	Orange	Pink
	Reference	3127251	3127252	3127253	3127254	3127255
Product picture						
Product description	Formatting easy	because armed with a copper metallic wire				
	Tips	removable, little, smooth, rounded				
	Use	ideal for all uses, even on irritated mucous membranes				
	Diametre	6.5 mm				
	Length	15 cm				
Characteristics and compositon	Made from PVC					
Packaging	<ul style="list-style-type: none"> - Store at room temperature in a dry and ventilated area - Plastic bag of 100 units 					
Shelf life	n.a					
Indications of use						
Legal information	Medical Device of class IIA CE 0051 Marking- Organismo de certificación: IMQ S.p.A Manufacturer: EURONDA					