



Product Designation	LATEX EXAMINATION GLOVES - LIGHTLY POWDERED														
Reference	XS -5/6 0116692	S - 6/7 0116656	M - 7/8 0116657	L-8/9 0116658	XL-9/10 0124592										
Product picture															
Product description	<ul style="list-style-type: none"> ▪ Powdered ambidextrous examination gloves, made of natural rubber latex. Non sterile, with smooth surface, with beaded cuff. 														
Characteristics and composition	<table border="1" data-bbox="579 931 1449 1055"> <thead> <tr> <th></th> <th>XS -5/6</th> <th>S - 6/7</th> <th>M - 7/8</th> <th>L-8/9</th> </tr> </thead> <tbody> <tr> <td>Palm Width (mm)</td> <td>70-79</td> <td>80-89</td> <td>90-99</td> <td>100-109</td> </tr> </tbody> </table> <ul style="list-style-type: none"> ▪ Protein content : max 2000 $\mu\text{g}/\text{dm}^2$ ▪ Powder level : max 10 mg/dm^2 ▪ Compliance standards : ASTM D3578-05, EN 455-1/2/3, EN 420, EN 388, EN 374-1/2/3 						XS -5/6	S - 6/7	M - 7/8	L-8/9	Palm Width (mm)	70-79	80-89	90-99	100-109
	XS -5/6	S - 6/7	M - 7/8	L-8/9											
Palm Width (mm)	70-79	80-89	90-99	100-109											
Packaging	Box containing 100 gloves (± 2 gloves). 10 boxes per carton.														
Shelf life	3 years from date of manufacture														
Indications of use															
Legal information	Medical Device of class I under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking														

Version 2.0
Updated 2014/04/22 by SMA

CFPM - 2 bis chemin du Loup 93290 Tremblay-en-France - FRANCE - Tel. +33 (0)1 45 91 31 82